


Thorlabs Deformable Mirror DMP40 Device Adapter

A brief user guide for Thorlabs deformable(DMP40) device adapter for Micro-Manager. This device adapter allows you to control the deformable mirror using device property browser through the tools tab.

Requirements

First of all you need to download Micro-Manager (https://micro-manager.org/wiki/Download%20Micro-Manager_Latest%20Release), and install it. we highly recommend download the nightly builds version. After installation, depend on the version of Micro-Manager, copy “mmgr_dal_Thorlab_Deformable_Mirror.dll” file to the Micro-Manager installation location file (C:\Program Files\Micro-Manager-1.4).

Follow instruction on: “https://micro-manager.org/wiki/Micro-Manager_Configuration_Guide” to configure your deformable mirror in Micro-Manager. You can add your mirror under ThorlabsDeformableMirror.



Usage

Thorlabs deformable mirror has properties that you can see on the image below.

1. Tip & Tilt: this property can change the 3 arm segments of mirror with amplitude and angle to modify the mirror tilt.
2. ApplyZernikes: this property only added for communication with REALM plugin.
3. Load Wavefront Correction: you can load segments voltage which you saved before for any purposes, for instance: flatness, tetrapod, saddle point and
4. Temperature: this property shows the mirror temperature.
5. Relax the Mirror: this property allows you to use relaxation feature offered by Thorlabs deformable mirror. Relax is a function that removes mechanical tensions in the mirror. For more information please read the Thorlabs deformable mirror manual.
6. Reset Segments: this property returns all 40 segments to 100V.
7. Reset Zernike: this property returns all Zernike coefficients to 0 and all 40 segments to 100V.
8. Save current position: this property save all segments voltages that you can use it later.
9. Mirror segments: these properties changes each segment voltage from 0V to 200V.
10. Zernike Coefficients: these properties changes each Zernike that you want. Moreover with this device adapter, it is designed so that if you change Zernike coefficients all segments will change depend on the each Zernike pattern.

Device Property Browser

- ☐ Show cameras
- ☐ Show shutters
- ☐ Show stages
- ☐ Show discrete changers
- ☒ Show other devices
- ☐ Show read-only properties
- Refresh!

Property	Value
ThorLab Deformable Mirror-Amplitude	0
ThorLab Deformable Mirror-Angle	0
ThorLab Deformable Mirror-ApplyZernikes	0
ThorLab Deformable Mirror-Load wavefront correction	DMDData/Init/WavefrontCorrection.txt
ThorLab Deformable Mirror-Mirror Temperature	29.1601
ThorLab Deformable Mirror-Relax the Mirror	OFF
ThorLab Deformable Mirror-Reset Segments	
ThorLab Deformable Mirror-Reset Zernikes	
ThorLab Deformable Mirror-Save current position [input filename]	DMDData/WavefrontCorrection_save.txt
ThorLab Deformable Mirror-Segment 1	100
ThorLab Deformable Mirror-Segment 10	100
ThorLab Deformable Mirror-Segment 11	100
ThorLab Deformable Mirror-Segment 12	100
ThorLab Deformable Mirror-Segment 13	100
ThorLab Deformable Mirror-Segment 14	100
ThorLab Deformable Mirror-Segment 15	100
ThorLab Deformable Mirror-Segment 16	100
ThorLab Deformable Mirror-Segment 17	100
ThorLab Deformable Mirror-Segment 18	100
ThorLab Deformable Mirror-Segment 19	100
ThorLab Deformable Mirror-Segment 2	100
ThorLab Deformable Mirror-Segment 20	100
ThorLab Deformable Mirror-Segment 21	100
ThorLab Deformable Mirror-Segment 22	100
ThorLab Deformable Mirror-Segment 23	100
ThorLab Deformable Mirror-Segment 24	100
ThorLab Deformable Mirror-Segment 25	100
ThorLab Deformable Mirror-Segment 26	100
ThorLab Deformable Mirror-Segment 27	100
ThorLab Deformable Mirror-Segment 28	100
ThorLab Deformable Mirror-Segment 29	100
ThorLab Deformable Mirror-Segment 3	100
ThorLab Deformable Mirror-Segment 30	100
ThorLab Deformable Mirror-Segment 31	100
ThorLab Deformable Mirror-Segment 32	100
ThorLab Deformable Mirror-Segment 33	100
ThorLab Deformable Mirror-Segment 34	100
ThorLab Deformable Mirror-Segment 35	100
ThorLab Deformable Mirror-Segment 36	100
ThorLab Deformable Mirror-Segment 37	100
ThorLab Deformable Mirror-Segment 38	100
ThorLab Deformable Mirror-Segment 39	100
ThorLab Deformable Mirror-Segment 4	100
ThorLab Deformable Mirror-Segment 40	100
ThorLab Deformable Mirror-Segment 5	100
ThorLab Deformable Mirror-Segment 6	100
ThorLab Deformable Mirror-Segment 7	100
ThorLab Deformable Mirror-Segment 8	100
ThorLab Deformable Mirror-Segment 9	100
ThorLab Deformable Mirror-Z2-2	0
ThorLab Deformable Mirror-Z20	0
ThorLab Deformable Mirror-Z22	0
ThorLab Deformable Mirror-Z3-1	0
ThorLab Deformable Mirror-Z3-3	0
ThorLab Deformable Mirror-Z3-3	0
ThorLab Deformable Mirror-Z31	0
ThorLab Deformable Mirror-Z33	0
ThorLab Deformable Mirror-Z4-2	0
ThorLab Deformable Mirror-Z4-4	0
ThorLab Deformable Mirror-Z40	0
ThorLab Deformable Mirror-Z42	0
ThorLab Deformable Mirror-Z44	0